


PICK1 Rabbit mAb

ABZ1810

Product Overview

Name	PICK1 Rabbit mAb
Accession(Primary)	Q9NRD5
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:20-1:50

Target information (Q9NRD5)

Gene Name PICK1**About Gene Symbol**

Western blot analysis of PICK1 in K562, HeLa lysates using PICK1 antibody.

