


MBNL1 Rabbit mAb

ABZ1811

Product Overview

Name	MBNL1 Rabbit mAb
Accession(Primary)	Q9NR56
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q9NR56)

Gene Name MBNL1**About Gene Symbol**

Western blot analysis of MBNL1 in A549, HL-60, U251 lysates using MBNL1 antibody.

