

DDX21 Rabbit mAb


ABZ1812

Product Overview


Name	DDX21 Rabbit mAb
Accession(Primary)	Q9NR30
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 87 kDa; Observed MW: 87 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:200-1:500; FC-1:50-1:100

Target information (Q9NR30)

Gene Name	DDX21
About Gene Symbol	



Immunocytochemistry analysis of DDX21 (green)
in 293T using DDX21 antibody, and DAPI (blue).



Western blot analysis of DDX21 in HeLa lysates
using DDX21 antibody.