

SMAC Rabbit mAb


ABZ1813

Product Overview


Name	SMAC Rabbit mAb
Accession(Primary)	Q9NR28
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 27 kDa; Observed MW: 21 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:20; ICC/IF-1:200-1:500; IP-1:10-1:100

Target information (Q9NR28)


Gene Name	DIABLO
About Gene Symbol	




Western blot analysis of Smac/Diablo in mouse heart, rat skin lysates using SMAC antibody.



Western blot analysis of Smac/Diablo in K562, HeLa lysates using Smac/Diablo antibody.



Immunocytochemistry analysis of Smac/Diablo(green) in HeLa using Smac/Diablo antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Smac/Diablo antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.