


TrpV5 Rabbit mAb

ABZ1818

Product Overview

Name	TrpV5 Rabbit mAb
Accession(Primary)	Q9NQA5
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 83 kDa; Observed MW: 83 kDa
Recommended Dilution	WB-1:2000-1:5000; IHC-P-1:100-1:200

Target information (Q9NQA5)

Gene Name TRPV5**About Gene Symbol**

Western blot analysis of TRPV5 in K562, C6, 3T3, HeLa lysates using TrpV5 antibody.

