

TMEM9B Rabbit mAb


ABZ1820

Product Overview


Name	TMEM9B Rabbit mAb
Accession(Primary)	Q9NQ34
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 23 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100

Target information (Q9NQ34)

Gene Name	TMEM9B
About Gene Symbol	TMEM9B (TMEM9 Domain Family Member B) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of TMEM9B in Hacat, PC-12 lysates using Transmembrane Protein 9B antibody.



Western blot analysis of TMEM9B in 3T3 lysates using TMEM9B antibody.