


## NOP10 Rabbit mAb

ABZ1822


### Product Overview

|                             |                                          |
|-----------------------------|------------------------------------------|
| <b>Name</b>                 | NOP10 Rabbit mAb                         |
| <b>Accession(Primary)</b>   | Q9NPE3                                   |
| <b>Host Species</b>         | Rabbit                                   |
| <b>Clonality</b>            | Monoclonal                               |
| <b>Isotype</b>              | IgG                                      |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                        |
| <b>Calculated MW</b>        | Calculated MW: 8 kDa; Observed MW: 8 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:200-1:500      |

### Target information (Q9NPE3)

**Gene Name** NOP10**About Gene Symbol**

Western blot analysis of NOP10 in HeLa, HL-60 lysates using NOP10 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using NOP10 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.