

PREB Rabbit mAb


ABZ1825

Product Overview


Name	PREB Rabbit mAb
Accession(Primary)	Q9HCU5
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 45 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200; IP-1:20-1:50

Target information (Q9HCU5)

Gene Name	PREB
About Gene Symbol	PREB (Prolactin Regulatory Element Binding) is an actin regulatory protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PREB in K562, C6, 3T3, HeLa lysates using PREB antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using PREB antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.