

LGR6 Rabbit mAb


ABZ1827

Product Overview


Name	LGR6 Rabbit mAb
Accession(Primary)	Q9HBX8
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 104 kDa; Observed MW: 104 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; ICC/IF-1:20-1:50; FC-1:10-1:100

Target information (Q9HBX8)


Gene Name	LGR6
About Gene Symbol	



Immunocytochemistry analysis of LGR6 (green) in HeLa using LGR6 antibody, and DAPI (blue).



Western blot analysis of GPCR LGR6 in K562, C6, 3T3, HeLa lysates using GPCR LGR6 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer tissue using GPCR LGR6 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.