

BARX1 Rabbit mAb

ABZ1828


Product Overview

Name	BARX1 Rabbit mAb
Accession(Primary)	Q9HBU1
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 27 kDa; Observed MW: 32 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q9HBU1)

Gene Name BARX1

About Gene Symbol



Western blot analysis of BarX1 in HeLa, CHO-K1, C6, Ramos lysates using BarX1 antibody.

