


TAPP1 Rabbit mAb

ABZ1829

Product Overview

Name	TAPP1 Rabbit mAb
Accession(Primary)	Q9HB21
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 46 kDa; Observed MW: 46 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:10-1:100

Target information (Q9HB21)

Gene Name PLEKHA1**About Gene Symbol**

Western blot analysis of TAPP1 in K562 lysates using PLEKHA1 antibody.

