

Nanog Rabbit mAb

ABZ1831

Product Overview

Name	Nanog Rabbit mAb
Accession(Primary)	Q9H9S0
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 35 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:50-1:100; IP-1:20-1:50

Target information (Q9H9S0)

Gene Name	NANOG
About Gene Symbol	