

PHD3 Rabbit mAb


ABZ1834

Product Overview


Name	PHD3 Rabbit mAb
Accession(Primary)	Q9H6Z9
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 27 kDa; Observed MW: 27 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:100-1:200; IP-1:50-1:100; IHC-1:50-1:100

Target information (Q9H6Z9)


Gene Name	EGLN3
About Gene Symbol	



Western blot analysis of PHD3 in HeLa lysates using PHD3 antibody.



Western blot analysis of PHD3 in rat lung and mouse lung lysates using PHD3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using PHD3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.