

EHD1 Rabbit mAb


ABZ1836

Product Overview


Name	EHD1 Rabbit mAb
Accession(Primary)	Q9H4M9
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 61 kDa; Observed MW: 61 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:1000; IP-1:10-1:100; ICC/IF-1:50-1:100

Target information (Q9H4M9)


Gene Name	EHD1
About Gene Symbol	



Immunocytochemistry analysis of EHD1 (green) in HeLa using EHD1 antibody, and DAPI (blue).



Western blot analysis of EHD1 in K562, rat Brain, C6, 3T3, HeLa lysates using EHD1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer tissue using EHD1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.