


TXNIP Rabbit mAb

ABZ1839

Product Overview

Name	TXNIP Rabbit mAb
Accession(Primary)	Q9H3M7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 44 kDa; Observed MW: 53 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:100-1:200

Target information (Q9H3M7)

Gene Name TXNIP**About Gene Symbol**

Western blot analysis of TXNIP in PC-12, NIH/3T3 lysates using TXNIP antibody.

Copyright© 2024 Zen-bioscience Limited.

