


FOXP1 Rabbit mAb

ABZ1840

Product Overview

Name	FOXP1 Rabbit mAb
Accession(Primary)	Q9H334
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 75 kDa; Observed MW: 82-90 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:500-1:1000

Target information (Q9H334)

Gene Name FOXP1**About Gene Symbol**

Western blot analysis of FOXP1 in Raw264.7, rat lung lysates using FOXP1 antibody.

