


OATP5A1 Rabbit mAb

ABZ1841


Product Overview

Name	OATP5A1 Rabbit mAb
Accession(Primary)	Q9H2Y9
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 92 kDa; Observed MW: 92 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100

Target information (Q9H2Y9)

Gene Name SLC05A1**About Gene Symbol**

Western blot analysis of SLC05A1 in T47D lysates using SLC05A1 antibody.



Western blot analysis of SLCO5A1 in K562, HeLa lysates using SLCO5A1 antibody.