

P2RY12 Rabbit mAb

ABZ1844

Product Overview

Name	P2RY12 Rabbit mAb
Accession(Primary)	Q9H244
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 39 kDa; Observed MW: 52 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:50-1:100

Target information (Q9H244)

Gene Name P2RY12

About Gene Symbol



