

Ctip1/BCL-11A Rabbit mAb


ABZ1846

Product Overview

Name	Ctip1/BCL-11A Rabbit mAb
Accession(Primary)	Q9H165
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 91 kDa; Observed MW: 100-145 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:20-1:100; ICC/IF-1:200-1:500; FC-1:100-1:200; IP-1:20-1:50

Target information (Q9H165)

Gene Name	BCL11A
About Gene Symbol	



Western blot analysis of Ctip1 in Ramos, Raji
lysates using Ctip1 antibody.