


ATG10 Rabbit mAb

ABZ1847

Product Overview

Name	ATG10 Rabbit mAb
Accession(Primary)	Q9H0Y0
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 25 kDa; Observed MW: 25 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q9H0Y0)

Gene Name ATG10**About Gene Symbol**

Western blot analysis of ATG10 in HeLa, A549, U251 lysates using ATG10 antibody.

