

NAT10 Rabbit mAb


ABZ1848

Product Overview


Name	NAT10 Rabbit mAb
Accession(Primary)	Q9H0A0
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 116 kDa; Observed MW: 116 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1:500-1:2000; FC-1:200-1:500; IP-1:50-1:100

Target information (Q9H0A0)

Gene Name	NAT10
About Gene Symbol	



Western blot analysis of NAT10 in Jurkat, C6, HeLa lysates using NAT10 antibody.



Immunocytochemistry analysis of NAT10(green) in HeLa using NAT10 antibody, and DAPI(blue)