


XPO4 Rabbit mAb

ABZ1853

Product Overview

Name	XPO4 Rabbit mAb
Accession(Primary)	Q9C0E2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 130 kDa; Observed MW: 130 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q9C0E2)

Gene Name XPO4**About Gene Symbol**

Western blot analysis of Exportin 4 in MCF-7, C6 lysates using XPO4 antibody.

