


Sprouty 4/Spry-4 Rabbit mAb

ABZ1856

Product Overview

Name	Sprouty 4/Spry-4 Rabbit mAb
Accession(Primary)	Q9C004
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 33 kDa; Observed MW: 33 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:10-1:100

Target information (Q9C004)

Gene Name SPRY4**About Gene Symbol**

Western blot analysis of Sprouty 4 in MCF-7, 3T3, C6 lysates using Sprouty 4 antibody.

