

BCL2L14/BCLG Rabbit mAb


ABZ1857

Product Overview


Name	BCL2L14/BCLG Rabbit mAb
Accession(Primary)	Q9BZR8
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 37 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:100-1:200; FC-1:50-1:100; IP-1:20-1:50

Target information (Q9BZR8)


Gene Name	BCL2L14
About Gene Symbol	



Western blot analysis of Bcl G/BCL2L14 in C6, 3T3, PC-12, Raw264.7 lysates using Bcl2 L14 antibody.



Western blot analysis of Bcl G/BCL2L14 in HeLa, A549, HL-60, U251, U87-MG lysates using Bcl G/BCL2L14 antibody.



Immunocytochemistry analysis of Bcl G/BCL2L14 (green) in hela using Bcl G/BCL2L14 antibody, and DAPI(blue)