


GTPBP4 Rabbit mAb

ABZ1858


Product Overview

Name	GTPBP4 Rabbit mAb
Accession(Primary)	Q9BZE4
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 74 kDa; Observed MW: 80 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500

Target information (Q9BZE4)

Gene Name GTPBP4**About Gene Symbol**

Western blot analysis of GTPBP4/NOG1 in A549, HL-60 lysates using GTPBP4 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using GTPBP4/NOG1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.