


ACE2 Rabbit mAb


ABZ1862

Product Overview


Name	ACE2 Rabbit mAb
Accession(Primary)	Q9BYF1
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 92 kDa; Observed MW: 120 kDa
Recommended Dilution	WB-1:500-1:1000; IHC:1:100-1:200

Target information (Q9BYF1)

Gene Name ACE2**About Gene Symbol**



Western blot analysis of ACE2 in Human kidney lysates using ACE2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human kidney. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.