

Angiopoietin-like 4/ANGPTL4 Rabbit mAb


ABZ1863

Product Overview

Name	Angiopoietin-like 4/ANGPTL4 Rabbit mAb
Accession(Primary)	Q9BY76
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 45 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:50-1:100; FC-1:200-1:1000

Target information (Q9BY76)

Gene Name	ANGPTL4
About Gene Symbol	



Western blot analysis of Angiopoietinlike 4 in K562, HeLa lysates using Angiopoietin Like 4 antibody.