

eIF2A Rabbit mAb


ABZ1864

Product Overview


Name	eIF2A Rabbit mAb
Accession(Primary)	Q9BY44
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 65 kDa; Observed MW: 65 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:50-1:100

Target information (Q9BY44)


Gene Name	EIF2A
About Gene Symbol	



Western blot analysis of eIF2A in 3T3, PC-12 lysates using eIF2A antibody.



Western blot analysis of eIF2A in LnCap, K562, U251MG lysates using eIF2A antibody.



Immunocytochemistry analysis of eIF2A (green) in MCF-7 using eIF2A antibody, and DAPI (blue)