


ITPA Rabbit mAb

ABZ1865


Product Overview

Name	ITPA Rabbit mAb
Accession(Primary)	Q9BY32
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 21 kDa; Observed MW: 21 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100

Target information (Q9BY32)

Gene Name ITPA**About Gene Symbol**

Western blot analysis of ITPA in HeLa lysates using ITPA antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Inosine triphosphate pyrophosphatase antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.