

PUMA Rabbit mAb


ABZ1867

Product Overview


Name	PUMA Rabbit mAb
Accession(Primary)	Q9BXH1
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 21 kDa; Observed MW: 21 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:200-1:500; FC-1:20-1:100

Target information (Q9BXH1)


Gene Name	BBC3
About Gene Symbol	




Immunocytochemistry analysis of PUMA (green) in 293T using PUMA antibody, and DAPI (blue).



Western blot analysis of PUMA in K562 lysates using PUMA antibody.



Western blot analysis of PUMA in Jurkat, THP-1, C6, 3T3, HeLa lysates using PUMA antibody



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using PUMA antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.