

## KIFC1 Rabbit mAb


ABZ1868

### Product Overview


|                             |                                                                                          |
|-----------------------------|------------------------------------------------------------------------------------------|
| <b>Name</b>                 | KIFC1 Rabbit mAb                                                                         |
| <b>Accession(Primary)</b>   | Q9BW19                                                                                   |
| <b>Host Species</b>         | Rabbit                                                                                   |
| <b>Clonality</b>            | Monoclonal                                                                               |
| <b>Isotype</b>              | IgG                                                                                      |
| <b>Species Reactivity</b>   | Human                                                                                    |
| <b>Calculated MW</b>        | Calculated MW: 74 kDa; Observed MW: 74 kDa                                               |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:1000-1:5000;<br>ICC/IF-1:50-1:100; FC-1:20-1:50; IP-1:10-1:100 |

### Target information (Q9BW19)


|                          |       |
|--------------------------|-------|
| <b>Gene Name</b>         | KIFC1 |
| <b>About Gene Symbol</b> |       |



Western blot analysis of KIFC1 in K562, 3T3, HeLa lysates using KIFC1 antibody.



Immunocytochemistry analysis of KIFC1(green) in HeLa using KIFC1 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using KIFC1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.