


PSS2 Rabbit mAb

ABZ1869

Product Overview

Name	PSS2 Rabbit mAb
Accession(Primary)	Q9BVG9
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 56 kDa; Observed MW: 56 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:20-1:100

Target information (Q9BVG9)

Gene Name PTDSS2**About Gene Symbol**

Western blot analysis of PTDSS2 in K562, HeLa lysates using PSS2 antibody.

