


## SPON2 Rabbit mAb

ABZ1870

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | SPON2 Rabbit mAb                           |
| <b>Accession(Primary)</b>   | Q9BUD6                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse                               |
| <b>Calculated MW</b>        | Calculated MW: 36 kDa; Observed MW: 36 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000                           |

### Target information (Q9BUD6)

**Gene Name** SPON2**About Gene Symbol**

Western blot analysis of SPON2 in K562, C6, 3T3, HeLa lysates using SPON2 antibody.

