

Transmembrane Protein 43 Rabbit mAb


ABZ1871

Product Overview


Name	Transmembrane Protein 43 Rabbit mAb
Accession(Primary)	Q9BTV4
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 45 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100

Target information (Q9BTV4)

Gene Name	TMEM43
About Gene Symbol	TMEM43 (Transmembrane Protein 43) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of TMEM43 in K562, C6 lysates using Transmembrane Protein 43 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using TMEM43 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.