

TRIM56 Rabbit mAb


ABZ1872

Product Overview


Name	TRIM56 Rabbit mAb
Accession(Primary)	Q9BRZ2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 81 kDa; Observed MW: 81 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:50-1:100; FC-1:20-1:50

Target information (Q9BRZ2)

Gene Name	TRIM56
About Gene Symbol	



Western blot analysis of TRIM56 in HeLa lysates using TRIM56 antibody.



Immunocytochemistry analysis of TRIM56(green) in HeLa using TRIM56 antibody, and DAPI(blue)