

ORAI3 Rabbit mAb

ABZ1873


Product Overview

Name	ORAI3 Rabbit mAb
Accession(Primary)	Q9BRQ5
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 31 kDa; Observed MW: 31 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:200-1:1000


Target information (Q9BRQ5)

Gene Name ORAI3

About Gene Symbol



Western blot analysis of Orai3 in HeLa, A549, HL-60, U251, U87-MG lysates using ORAI3 antibody.



Immunohistochemistry analysis of paraffin-embedded liver cancer using Orai3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.