


PCGF6 Rabbit mAb

ABZ1874

Product Overview

Name	PCGF6 Rabbit mAb
Accession(Primary)	Q99NA9
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 40 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20

Target information (Q99NA9)

Gene Name Pcgf6**About Gene Symbol**

Western blot analysis of PCGF6 in PC-12 lysates
using PCGF6 antibody.

