

GDF15 Rabbit mAb

ABZ1875


Product Overview

Name	GDF15 Rabbit mAb
Accession(Primary)	Q99988
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q99988)

Gene Name GDF15

About Gene Symbol



Western blot analysis of GDF15 in HeLa lysates using GDF15 antibody.

