

BAG1 Rabbit mAb


ABZ1876

Product Overview


Name	BAG1 Rabbit mAb
Accession(Primary)	Q99933
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 39 kDa; Observed MW: 52,46,33 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200

Target information (Q99933)

Gene Name	BAG1
About Gene Symbol	BAG1 (BAG Cochaperone 1) is a molecular chaperone involved in apoptosis. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Bag1 in HeLa, MCF-7, U251MG, K562 lysates using BAG1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using Bag1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.