

TSG101 Rabbit mAb


ABZ1878

Product Overview


Name	TSG101 Rabbit mAb
Accession(Primary)	Q99816
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 44 kDa; Observed MW: 44 kDa
Recommended Dilution	WB-1:2000-1:20000; IHC-1:100-1:1000; IF-1:200-1:500; FC-1:20-1:50

Target information (Q99816)


Gene Name	TSG101
About Gene Symbol	



Immunocytochemistry analysis of TSG101 (green) in HT-1080 using TSG101 antibody, and DAPI (blue).



Western blot analysis of TSG101 in K562, rat Brain, C6, 3T3, HeLa lysates using TSG101 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil tissue using TSG101 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.