

Smad5 Rabbit mAb


ABZ1881

Product Overview


Name	Smad5 Rabbit mAb
Accession(Primary)	Q99717
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 52 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:20-1:50; ICC/IF-1:50-1:100; FC-1:20-1:100

Target information (Q99717)

Gene Name	SMAD5
About Gene Symbol	SMAD5 (SMAD Family Member 5) is a transcription factor involved in TGF- β signaling. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic ta



Western blot analysis of Smad5 in C6, HeLa lysates using Smad5 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using SMAD5 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.