

ASK1 Rabbit mAb


ABZ1882

Product Overview

Name	ASK1 Rabbit mAb
Accession(Primary)	Q99683
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 155 kDa; Observed MW: 155 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:20-1:50; ICC/IF-1:100-1:200; FC-1:50-1:100

Target information (Q99683)

Gene Name	MAP3K5
About Gene Symbol	



Western blot analysis of ASK1 in HeLa lysates
using ASK1 antibody.