

LGMN Rabbit mAb


ABZ1885

Product Overview

Name	LGMN Rabbit mAb
Accession(Primary)	Q99538
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 49 kDa; Observed MW: 56,46,37 kDa
Recommended Dilution	WB-1:2000-1:5000

Target information (Q99538)

Gene Name	LGMN
About Gene Symbol	LGMN (Legumain) is a protease enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Legumain in Jurkat, HeLa lysates using Legumain antibody.