

PARK7/DJ1 Rabbit mAb


ABZ1887

Product Overview


Name	PARK7/DJ1 Rabbit mAb
Accession(Primary)	Q99497
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 20 kDa; Observed MW: 20 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:50-1:500; FC-1:20-1:50; IP-1:20

Target information (Q99497)


Gene Name	PARK7
About Gene Symbol	PARK7 (Parkinsonism Associated Deglycase) is a metabolic enzyme involved in cellular metabolism. It is an important research target in cancer research and neurodegenerative disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PARK7/DJ1 in K562, rat Brain, C6, 3T3, HeLa lysates using PARK7/DJ1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using PARK7/DJ1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of PARK7/DJ1
(green) in Hela using PARK7/DJ1 antibody, and
DAPI(blue)