


NOTCH4 Rabbit mAb

ABZ1889

Product Overview

Name	NOTCH4 Rabbit mAb
Accession(Primary)	Q99466
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 210 kDa; Observed MW: 36 kDa
Recommended Dilution	WB-1:2000-1:5000

Target information (Q99466)

Gene Name NOTCH4**About Gene Symbol**

Western blot analysis of NOTCH4 in HeLa, U87-MG lysates using Notch4 antibody.

