

MUC4 Rabbit mAb


ABZ1891

Product Overview

Name	MUC4 Rabbit mAb
Accession(Primary)	Q99102
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 542 kDa; Observed MW: 120 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q99102)

Gene Name	MUC4
About Gene Symbol	MUC4 (Mucin 4, Cell Surface Associated) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of MUC4 in PC-3 lysates using Mucin 4 antibody.