

HBEGF Rabbit mAb

ABZ1892


Product Overview

Name	HBEGF Rabbit mAb
Accession(Primary)	Q99075
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 23 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:20-1:50

Target information (Q99075)

Gene Name HBEGF

About Gene Symbol



Western blot analysis of HBEGF in K562, C6, 3T3, HeLa lysates using HBEGF antibody.

