

UHRF1 Rabbit mAb


ABZ1893

Product Overview

Name	UHRF1 Rabbit mAb
Accession(Primary)	Q96T88
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 90 kDa; Observed MW: 97 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:20-1:100; ICC/IF-1:200-1:500; FC-1:50-1:100; IP-1:20-1:50

Target information (Q96T88)

Gene Name	UHRF1
About Gene Symbol	



Western blot analysis of UHRF1 in K562, HeLa lysates using UHRF1 antibody.