

## NODAL Rabbit mAb

ABZ1895


### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | NODAL Rabbit mAb                           |
| <b>Accession(Primary)</b>   | Q96S42                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Rat                                 |
| <b>Calculated MW</b>        | Calculated MW: 40 kDa; Observed MW: 40 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; FC-1:100-1:500           |

### Target information (Q96S42)

**Gene Name** NODAL

**About Gene Symbol**



Western blot analysis of Nodal in C6, 3T3, HeLa lysates using Nodal antibody.

