

TRIB3 Rabbit mAb


ABZ1896

Product Overview

Name	TRIB3 Rabbit mAb
Accession(Primary)	Q96RU7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 40 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:10

Target information (Q96RU7)

Gene Name	TRIB3
About Gene Symbol	TRIB3 (Tribbles Pseudokinase 3) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of TRIB3 in K562, HeLa
lysates using TRIB3 antibody.